



STANDING SLING with 2 belts

A unique standing sling, which supports user when lifting from a stand aid.

USAGE

- For users with good upper body control and some lower limb function.
- Used for transfers from a sitting position.
- Used for transfers in stand aids.

FEATURES

- The sling is for transferring user from a wheelchair/chair.
- The sling provides good support around the back, giving the user a sense of security, when lifted and when standing.
- The sling with two belts can be adapted to users with a big belly, also the belt can be tightened over and under the belly.
- Standing sling has anti-slip on the inner side of the sling. This ensures that the sling does not slip up during lifting.
- The sling has padding under the arms.

SPECIFICATIONS

SWL: 250 kg

Material:
Polyester

Sizes: XS-S-M-L-XL-XXL

ITEM NUMBERS:

Standing sling with two belts, polyester:

	Item no.
XL	009-09563
S	009-09564
M	009-09565
L	009-09566
XL	009-09567
XXL	009-09568




SLING SIZE GUIDE

The sling has a coloured strap on upper back with size and midline mark, which ensures correct placement of the sling.

Colour code is presented on the size chart.

The image below shows a yellow colour code, which means that the sling is in size M (Medium):



	XS	S	M	L	XL	XXL
Belly dimension	62 – 80 cm	66 – 90 cm	80 – 100 cm	100 – 126 cm	110 – 132 cm	120 – 140 cm



Please consult the technical manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. Winncare Nordic reserves the right to changes without notice. Product photos may include accessories. Please consult your dealer for further information.

FOR MORE INFORMATION VISIT

winncare.com

E-mail: info@winncare.dk

Phone: +45 70273720

