



COMFORT DISPOSABLE SLING

Non-woven disposable full sling. Available with or without head support.

USAGE

- A disposable sling is perfect for situations where a high hygienic standard is needed, e.g. if there is a risk of infection.
- For users with increased disability or other special needs.
- For transfers from a sitting or lying position.
- For transfers in mobile and stationary hoists.

FEATURES

- The sling provides good support of the body and stabilizes the pelvis, supporting the upper body and giving the user a feeling of comfort and safety.
- The leg straps have a large padded surface, which reduces the pressure on the entire femur.
- The sling ensures a good upright posture that facilitates the positioning of the user.
- The sling is made from non-woven material.
- Available with or without head support.

SPECIFICATIONS

SWL: 275 kg

Material: PP, Non-woven

Sizes: S-M-L-XL

Available with or without head support.

ITEM NUMBERS

Disposable Comfort Sling without head support (10 pcs.)

Item numbers	
S	008-04104
M	008-04105
L	008-04106
XL	008-04107

Disposable Comfort Sling with head support (10 pcs.)

Item numbers	
S	008-04114
M	008-04115
L	008-04116
XL	009-04117



PRODUCT HIGHLIGHTS

- Split leg straps, which do not cross, due to an extra closing strap.
- Small seat opening.
- Washing symbol with an indicator which warns the carer if the sling is damp/wet and therefore needs to be destroyed and replaced:



- A coloured strap on upper back with size and midline mark, which ensure correct placement of the sling. Colour code is presented on the size chart. The image below shows a yellow colour code, which means that the sling is in size M (Medium):



	S	M	L	XL
Y Hip measurement	42 cm	48 cm	52 cm	68 cm
X Without head support Measurement from head to lower back.	55 cm	66 cm	75 cm	75 cm
Z With head support: Measurement from neck to lower back.	80 cm	91 cm	100 cm	100 cm



Please consult the technical manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. WinnCare Nordic reserves the right to changes without notice. Product photos may include accessories. Please consult your dealer for further information.



FOR MORE INFORMATION, VISIT

winnCare.com
E-mail: info@winnCare.dk
Phone: +45 70273720

WIN  CARE