



UNIVERSAL HAMMOCK

Unique sling with entire body support and good comfort

USAGE

- For less active users, who are in seated position for longer periods.
- For transfers from a lying position.
- For transfers in mobile and stationary hoists.

FEATURES

- The sling is transferring user from lying position or from the bed or bath.
- The sling provides good support of the body, which gives the user a feeling of safety during transfers.
- The sling ensures an upright sitting position that contributes to a correct positioning of the user.
- The sling is available in polyester mesh, which is perfect for bathing.
- The sling has unsplit leg straps.

SPECIFICATIONS

SWL: 275 kg

Material:

Polyester mesh.

Sizes: XS-S-M-L-XL-XXL

ITEM NUMBERS:

Universal Hammock, polyester mesh:

| | Item no. |
|------------|-----------|
| XS | 009-09673 |
| S | 009-09674 |
| M | 009-09675 |
| L | 009-09676 |
| XL | 009-09677 |
| XXL | 009-09678 |





SLING SIZE GUIDE

The sling has a coloured strap on upper back with size and midline mark, which ensures correct placement of the sling.

Colour code is presented on the size chart.

The image below shows a yellow colour code, which means that the sling is in size M (Medium):



| | XS | S | M | L | XL | XXL |
|--|-------|-------|-------|-------|-------|-------|
| Y (hip measurement) | 37 cm | 42 cm | 48 cm | 52 cm | 68 cm | 78 cm |
| X (measurement from neck to lower back) | 50 cm | 55 cm | 66 cm | 75 cm | 75 cm | 75 cm |



Please consult the technical manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. WinnCare Nordic reserves the right to changes without notice. Product photos may include accessories. Please consult your dealer for further information.

FOR MORE INFORMATION VISIT

winnCare.com

E-mail: info@winnCare.dk

Phone: +45 70273720

